

## Standard Features:

### Construction

- 6 storage drawers in center
- Slide-out keyboard tray
- Foot pedal with water on/off toggle
- Choice of Silestone® work surface and custom paint colors
- Instrument swivel mount bar
- Instrument air & water convenience outlets
- Individual handpiece pressure and coolant adjustments
- Rail system storage with 2 storage bins
- Solids trap
- Internal power strip on back for accessory items
- Aircraft grade aluminum chassis construction
- Baked on powder coat finish for durability
- Hospital grade plug and listed wiring and devices
- Rear storage compartment

### Handpieces

- 1 high speed air handpiece connection
- 1 low speed air handpiece connection
- Titanium straight HVE
- Titanium saliva ejector
- Air/water syringe

### Water System

- Single 1-liter self-contained water system
- DentaPure® water system disinfection – a single cartridge provides water treatment for one full year

### Casters

- Commercial grade 3" dual-wheel casters

### Utilities

- Master air supply regulator and internal air filter
- Umbilical with compressed air line, 1 main power cord, and 1 dedicated accessory power cord

