

MOMENTUM SADDLE™

MODEL: 90-3508



ASI's Momentum Saddle™ Model 90-3508

MAINTAIN YOUR DAILY MOMENTUM

ASI's Momentum stools are ergonomically designed combining comfort, optimal support and freedom of movement to help keep you going all day.

SIT OR STAND

The Momentum Saddle™ is designed to be used in sitting and semi-standing positions.

The saddle shape of the seat allows for almost full leg extension, making it easier to sit higher without putting pressure on the thighs.



The Leader In Advanced Systems Integration

Visit Us Online: www.asiDental.com

Email: sales@asidental.com | Toll Free: 844.880.3636 | Local: 303.407.6073

FEATURES & SPECIFICATIONS

ERGONOMIC SUPPORT

A good ergonomic stool needs to have a lumbar support to prevent slouching. The small lumbar support in the Momentum Saddle™ prevents flattening of the lower back, allowing freedom of movement and stretching of the spine. It is proven that we lose at least 30% of the curvature of the spine after 5-10 minutes of sitting. Without support, more pressure will occur within the discs increasing the chance of spinal injury and lower back problems.

STANDARD FEATURES:

Adjustable backrest for lumbar support. Slip-resistant Valencia® upholstery is coated with a PermaBlock™ protective finish for easy cleaning. Seat tilt and height adjustment for maximum operator comfort.



ASI's Momentum Saddle™
Model 90-3508

AVAILABLE OPTIONS

Foot Ring	For operators who need to sit higher
Tall Cylinder	For operators over 6' tall
Custom Colors	To view color options, visit https://vyvafabrics.com/vinyl/valencia/

SPECIFICATIONS

Cylinder Height	Standard
Standard Solid Color Options	Taupe, Bronze, Mushroom (Brown), Praline (Chocolate), Titan (Steel Blue), Platin (Gray), Graphite (Dark Gray), Schwarz (Black)
Standard Two-Tone Color Options	Schwarz with Platin Accent Edge, Taupe with Bronze Accent Edge