

# MOMENTUM MAGNIFICATION™

MODEL: 90-3512



*ASI's Momentum Magnification™*

## MAINTAIN YOUR DAILY MOMENTUM

ASI's Momentum stools are ergonomically designed for comfort, optimal support and freedom of movement to help keep you going all day.

## CUSTOM FIT & BALANCE

ASI's Momentum Magnification™ chair features an adjustable front seat angle that works in conjunction with the convex shaped back rest for a custom fit.

Achieve your preferred balance of support from the thighs to the back, without rounding or overloading the spine.



Visit Us Online: [www.asidental.com](http://www.asidental.com)

Email: [sales@asidental.com](mailto:sales@asidental.com) | Toll Free: 844.880.3636 | Local: 303.407.6073

8811 American Way, Suite 130 | Englewood, CO 80112



# MOMENTUM MAGNIFICATION™



## FEATURES & SPECIFICATIONS

### ERGONOMIC SUPPORT

A common tendency is to sit forward on the edge of your seat. This can cut off circulation to the legs and leave the back unsupported, resulting in a backward rounding of the spine. The hydraulically controlled front seat angle can be positioned downward to provide support. Used in conjunction with the seat tilt, the operator can sit higher. This position increases the hip angle to a more open 110°, reducing pressure on the spine. This allows the operator to remain fully seated at the back of the chair, receiving optimal back support while maintaining excellent proximity to patients.

The back rest is shaped to support the lumbar curve. Adjusted in conjunction with the seat angle, it provides support for long procedures and prevents slouching from muscle fatigue. Rubber bushings allow the backrest to swivel vertically, which aids in stretching while supporting the back fully.

### STANDARD FEATURES

Adjustable backrest for lumbar support. Slip-resistant Valencia® upholstery is coated with a PermaBlock™ protective finish for easy cleaning. Seat tilt and height adjustment for maximum operator comfort.

### AVAILABLE OPTIONS

<b>Foot Ring</b>	For operators who need to sit higher
<b>Tall Cylinder</b>	For operators over 6' tall
<b>Custom Colors</b>	To view color options, visit <a href="http://yvafabrics.com/vinyl/valencia">yvafabrics.com/vinyl/valencia</a>

### SPECIFICATIONS

<b>Cylinder Height</b>	Standard
<b>Standard Solid Color Options</b>	Taupe, Bronze, Mushroom (Brown), Praline (Chocolate), Titan (Steel Blue), Platin (Gray), Graphite (Dark Gray), Schwarz (Black)
<b>Standard Solid Color Options</b>	Schwarz with Platin Accent Edge, Taupe with Bronze Accent Edge



Visit Us Online: [www.asidental.com](http://www.asidental.com)

Email: [sales@asidental.com](mailto:sales@asidental.com) | Toll Free: 844.880.3636 | Local: 303.407.6073

8811 American Way, Suite 130 | Englewood, CO 80112

