

WHEN
PERFORMANCE
MATTERS

ADVANCED
ENDODONTIC
SYSTEMS[®]

FLEXIBLE TREATMENT ROOM DESIGN
ADAPTABLE TO THE WAY YOU WORK
ENHANCE PATIENT EXPERIENCE

PRECISION DESIGN

FOR THE WAY YOU PRACTICE

BUILD A MINIMALIST OFFICE WITH UNCONFINED TREATMENT ROOMS

ASI's space-saving designs create different ways to setup treatment suites that do not exist with traditional side and rear cabinetry. Our systems adapt to a number of room layouts, sizes, and dual or single-entry scenarios. Their mobility solves ambidextrous challenges.

Our total solution package is a flexible, cart-based delivery for both Doctor and Assistant. The modular platform creates freedom in design, providing superior ergonomic function and the ability to update with newly emerging technology.

Plug-n-play functionality means adding equipment only as necessary. Grow your practice at your own pace while avoiding additional construction & equipment costs.

REDEFINE & IMPROVE PATIENT ACCEPTANCE

Give your patients the unexpected. Allow them to feel welcome, safe, and confidently reassured by eliminating the overwhelming clutter of chair mounted arms, instrument control boxes, and tangled power and computer cables. Present an image of specialization that comes with the professional and unique integration that an ASI Advanced Endodontic System® provides. Patients want the assurance and peace of mind that they are receiving the most advanced treatment possible. A neatly organized, streamlined system with an advanced medical technology appearance reinforces their confidence and their experience with you.



Anthracite
Metallic



Newport

(samples pictured on system)



PRECISION EQUIPMENT THAT DELIVERS

CUSTOMIZABLE INTEGRATION

A comprehensive design that provides everything you need right at your fingertips. The ready to go delivery system includes primary instrumentation of pneumatic handpiece connections, air/water, and endo air syringe. ASI's adaptable modular instrument panels and holder bar can be integrated with endodontic instruments of your choice in the order you desire. The systems are built for the way you practice, with only the instrumentation you require, allowing you to have as little or as many instruments you need with a single foot control that streamlines operation. Select ASI under floor foot control option for ease of mobility.


Simplified field upgradeability allows instruments to be replaced in the future with newer models right in your office. Improve your practice efficiency and ergonomics by eliminating the clutter and confusion of table top instruments with separate power cords and pedals. Have what you need within ease of reach including ASI unique endodontic instruments.

- Endodontic Piezo Ultrasonic with saline irrigation option
- Electric Endo Motors and Surgical Motors
- Micro Irrigation Instrumentation
- Cordless Instrument Mounting with integrated power station
- Illumination Instruments
- IT integration capable with data cabling and monitor mounting

HIGHER QUALITY AND ADVANCED DESIGN FEATURES

- High-grade ASI MetalCraft™ aluminum architecture for strength and lifetime use
- Powder coat finish and engineered non-porous work tops that are resilient to disinfectants and sodium hypochlorite
- ASI White Silk Handpiece Tubing that is soft, non-sticking and light to use
- Stainless steel and titanium components and valves
- ASI TruFit™ Articulating Handpiece holders for the perfect fit for each instrument with personalized placement
- Integral water line disinfection system
- Made in USA with Safety Standard Certifications

Customize Your System with
different paint and work top options



(samples pictured on system)

go to ASIDental.com/Electric-Motors-Instruments to browse our selection of Electric Motors and Instruments to integrate into your ASI Endodontic System.

WE'VE EARNED OUR QUALITY REPUTATION...

DESIGNER SERIES DELIVERY SYSTEMS

Versatile and full-featured. The ASI Designer Series combines modern esthetic styles with enhanced integration capabilities. Style features include high-visual angulation for control displays, the ability to choose from custom color selections, and engineered stone worktops to make the perfect complement to your treatment room. Increased in-cart organization of additional instruments and Information Technology (IT) components for the widest range of electronics and treatment capabilities including surgical motors and external irrigation. Select from a complete line of models for the right need. View color and stone worktop selections at ASIDental.com.

iTech Model: Seamless concealed computer cabling and equipment (90-2134)

Drawer Model: Includes five storage drawers for supplies (90-2133)

Compressor Model: Side delivery system – no utilities – just plug it in (90-2115)

Dual Compressor & Vacuum Model: Ideal for rear treatment delivery – no utilities – just plug it in (90-2125)

DESIGNER SUPPORT WORKSTATION

Advanced delivery for the dental assisting support team with quality titanium suction instruments that are impervious to bleach and harsh chemicals. Comprehensive design creates a flexible workstation that combines suction, optional nitrous oxide, and irrigation delivery with built-in storage, an ample engineered stone worktop, and an array of customizable supply organizers including optional instrument holders. IT integration provides sliding monitor mount, internal CPU compartment, integrated cabling with side mounted USB connections, and slide out keyboard tray. All mounted on a convenient mobile platform that easily moves within the treatment area so that the dental assistant team member can optimize for their most efficient and ergonomic use.

Ambidextrous Rear Delivery: Carriage mounted sliding instrument holder arm can be positioned anywhere along the radius track for preferred placement. (90-1054)



... LET US HELP ESTABLISH YOURS



CLASSIC SERIES DELIVERY SYSTEMS

Shorter with a more compact footprint allows the system to fit well in smaller treatment rooms and can be readily moved to either side of the treatment room for ambidextrous requirements. Compact but capable allows full integration of endodontic instrumentation and IT integration.

Door Model: Enclosed slide out water supply (90-2143E)

Drawer Model: Includes storage drawers for supplies (90-2140E)

ReFit™ Model: Reduced height to fit under existing dental countertops and to make an ideal retrofit to replace older dental equipment with a fresh, modern design (90-2143ES)

CLASSIC SUPPORT WORKSTATION

Advanced delivery for the dental assisting support team with quality titanium suction instruments that are impervious to bleach and harsh chemicals.

Freedom Model: The right size to fit most back wall areas while giving curved instrument delivery and to allow the cart to move closer to the treatment area. Features a slide out tray and supply organization (90-1044)

Mobile Assistant Cart Model: Small compact storage work cart with swivel arm mounted suction instruments and storage drawers (90-1021A)


Ambidextrous Assistant Cart Model: Same as 90-1021A Model above but with ability to relocate assistants' instruments arms to either side of the cart chassis (90-1021X)



go to ASIDental.com to
see our full product line
to meet your needs!

UNIQUE DESIGN FEATURES NOT FOUND ELSEWHERE

DON'T SETTLE FOR LESS. **DESIGN THE SYSTEM THAT MEETS YOUR NEEDS.**



LEADER IN ENDODONTIC INNOVATION

In the mid-1990s, ASI began working with Dr. Gary Carr in developing advanced delivery systems designed for the unique needs of the endodontic specialty. Since then, ASI has been committed to devoting time and resources to the ongoing development of new and improved instruments and delivery systems that optimize the workflow for Endodontists. We highly value and utilize customer input combined with our engineering expertise to create innovations, and will continue to do so for our customers. We specifically manufacture our products here in the U.S. not only so that we can readily control the quality of manufacturing, but also easily modify and improve designs as desired. Further, we endeavor to make products of superb quality by using high-grade materials including aluminum architecture, stainless steel and titanium pneumatic valves, durable finishes, and durable hardware. You need to experience an ASI Delivery System to compare the look, feel and quality.