

When Precision Matters

DESIGNED FOR THE ORAL SURGEON
AND IMPLANT SPECIALIST




Proudly Made
in the USA!

844.880.3636 | ASIDental.com
8811 American Way, Suite 130 | Englewood, CO 80112

 **ASI**
Dental Specialties

Advanced Dental Systems®

The Glider™

INTEGRATED DENTAL SURGICAL TABLE


Model 90-2148



Proudly Made
in the USA!

An integrated instrument solution for efficient implant and surgical procedures. The Glider™ is easily positioned over standard and surgical style patient chairs and provides sufficient workspace to lay out sterile instruments while integrating surgical motors, handpieces, and suction instruments for ease of access by the dental team.

The Glider is mounted on stable, heavy-duty casters with a programmable electric height adjustment. This allows for a full range of operator capability for either sit-down or stand-up treatment with mobile positioning for desired placement.



STRUCTURAL FEATURES


- Smooth, easy-to-clean, non-porous solid surface worktop provides sufficient workspace for placing hand instruments and other devices
- Telescoping height-adjustable with programmable electric adjustment
- High-grade dual wheel casters for strength and ease of mobility
- Aluminum MetalCraft design for esthetics, durable lifetime quality and lighter weight
- Pull out center drawer with locking option
- Rear air and water quick connects and electrical outlets for table top instruments
- Rear-mounted holder for disposable suction canisters or solids trap filter for connection to central suction

INSTRUMENT DELIVERY FEATURES

- Customizable instrument delivery with operation from only one foot control
- Standard pneumatic air handpiece connections and dual air/water syringes for doctor and assistant
- Titanium suction handpieces with autoclavable tubing option
- Articulating handpiece holder bar for preferred placement
- Internal water bottle system for standard irrigation with surgical saline pump option

INSTRUMENT & ACCESSORY OPTIONS

- Integrated surgical motors with full range of handpiece attachments for surgical to implant
- Integrated piezo ultrasonic options
- Integrated curing light option
- Side-mounted dovetail bars on each side. Holds accessory items including saline hooks, sharps container, removable sterilizable handpiece holder, wipes bracket, and waste bins
- Side-mounted or rear mounted glove dispenser option
- Rear rail system for storage trays, glove box holders, or flip top supply bins
- Post-mounted monitor mounts with either VESA mount or shelf



SCAN NOW FOR FULL
SPECIFICATIONS & PRICING

844.880.3636 | ASIDental.com

Implant to Surgical

EQUIPMENT SOLUTIONS


IMPLANT-SURGICAL PROCEDURE CART

High-quality modular procedure cart provides treatment room flexibility and is designed to complement your dental treatment equipment set up.



INTEGRATED DENTAL SYSTEM

Fully integrated delivery system provides convenient access to desired instruments with one foot control and without the clutter of separate control boxes and cords.



SCAN NOW
FOR FULL
SPECIFICATIONS
& PRICING



Proudly Made
in the USA!

844.880.3636 | ASIDental.com
8811 American Way, Suite 130 | Englewood, CO 80112